drive

 DeVilbiss®
HEALTHCARE

CONTENTS.....1

FINGERTIP PULSE OXIMETER.....2

FINGER-PULSOXIMETER.....20

PULSIOXÍMETRO.....41

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PULSOSSIMETRO DA DITO.....82

English

OPERATOR'S MANUAL

Rev B

FINGERTIP PULSE OXIMETER

CE
0123

General Description

Oxygen saturation (SpO_2) is a percentage of oxygen contained in the blood. 'S' indicates saturation, 'p' indicates the pulse and O_2 is the symbol for oxygen. Blood cells contain a chemical called hemoglobin (Hb). When hemoglobin is combined with oxygen it forms a second chemical called oxyhemoglobin (HbO_2). As blood travels around the body the oxyhemoglobin releases the oxygen where required and then turns back into hemoglobin. Oxygen levels in the blood can fluctuate, especially during activity. Generally normal oxygen levels of a

person is 95% or greater. Less than 95% may indicate health issues which would require clinical investigation. Oxygen deficiency in the blood is referred to as hypoxia. There are many respiratory conditions that can result in a patient becoming hypoxic; COPD (Chronic Obstructive Pulmonary Disease), OSA (Obstructive Sleep Apnea) and chronic bronchitis to name a few. Oxygen saturation plays an important part for clinical professionals to be proactive in patient diagnosis.

Perfusion Index (PI) is related to the strength of the patients pulse at the site of measurement. PI is measured as a percentage (%) and the optimal value is 20% indicating a very strong pulse.

The drive DeVilbiss Healthcare HbO-3000 fingertip pulse oximeter is compact in size, user friendly and features an OLED screen display. SpO₂ diagnosis

requires the patient to simply wear the device on their finger. The oxygen saturation and pulse rate values will be displayed on the device screen. The drive DeVilbiss Healthcare HbO-3000 has been designed and tested to provide repetitive, accurate screening for SpO₂ patient diagnosis.

Measurement Principle

The principle measurement of an oximeter uses a mathematical formula known as the 'Lambert Beer Law' according to spectrum absorption characteristics of reductive hemoglobin (RHb) and Oxyhemoglobin (HbO₂) in glow and near-infrared zones. Photoelectric Oxyhemoglobin Inspection Technology is adopted in accordance with capacity pulse scanning and recording technology, so that two beams of lights at different wavelengths (660nm

glow and 905nm near infrared light) can be applied to the fingertip via a clamp finger-type sensor. A measured signal obtained by a photosensitive element, and processed via the device electronics and microprocessor will be displayed on the oximeter's screen display.

Diagram of Operation Principle

- 1.Red and Infrared-ray Emission Tube
- 2.Red and Infrared-ray Receipt Tube

Product Operation Scope

The HbO-3000 fingertip pulse oximeter is a non-invasive device designed for oxygen saturation (SpO_2) and pulse rate detection of adult and pediatric patients in the home and clinical environments.

It has not been designed for continuous monitoring.

Precautions for use

1. Carefully read this user manual before using the HbO-3000 fingertip pulse oximeter.
2. Do not use the HbO-3000 fingertip pulse oximeter in an MRI or CT environment.
3. A pulse oximeter is intended only as a supplementary accessory in patient assessment. It must be used in conjunction with other methods of clinical diagnosis.
4. Frequently check application site to determine correct positioning of the sensor.
5. Check patients circulation before use and refrain from using if the skin becomes sensitive.
6. Prolonged use of the device or the

patient's condition may require changing the site of the sensor periodically.

7. Autoclaving, ethylene oxide sterilizing, or immersing the HbO-3000 fingertip pulse oximeter in liquids will damage the device.
8. The following conditions could cause inaccurate readings:
 - High levels of dysfunctional haemoglobins such as carboxyhaemoglobin or methemoglobin.
 - Intravascular dyes such as indocyanine green or methylene blue.
 - High ambient light.
 - Excessive patient movement.
 - High-frequency electrosurgical interference.
 - Venous insufficiency.
 - Application of a blood pressure cuff, arterial catheter, or intravascular line.
 - Hypotension, severe vasoconstriction,


severe anemia, or hypothermia.

- Cardiac arrest or shock.
 - Fingernail polish or cosmetic fingernails.
- Please comply with local authority regulations when disposing/recycling of the device and batteries.

Product features

1. The HbO-3000 Fingertip pulse oximeter is light weight, compact and easy to use.
2. Low power consumption only requires 2 x AAA batteries. (20 hours).
3. A low voltage warning will be displayed on the device screen when batteries are low.
4. The product will automatically power off after 8 seconds if no signal is detected.
5. No maintenance required.

Operation Instructions

The HbO-3000 features six different directional screen settings including a pulse rate bar, plethysmography waveform and a PI indicator. In addition, the HbO-3000 has 5 different brightness levels for user preference.

Parameter Settings

Press the power button (>0.5s), the oximeter will enter into the parameter setting screen. There are two parameter sub menu screens:

Devices with sound/alarm function

When the “*” signal is shown on the “Sounds Setup”, press the button(>0.5s)

and enter into the sounds setting menu, pressing the button in turn will select the item to be changed, Press and hold the button to change the parameter. Select “+” or “-” to increase/decrease the number values.

When the “*” signal is shown on the “Alm Setup”, press the button (>0.5s) and enter into the alarm setting menu. Press the button to set the on/off for the alarm and beep.


When the “*” signal is shown on the “Restore”, press the button (>0.5s) to default to factory settings.

The ‘Demo’ feature is for demonstration purposes only.

Settings	
Alm Setup *	
Alm	on
Beep	off
Demo	off
Brightness	4
Restore	OK
Exit *	

Settings	
Sounds Setup *	
SpO ₂ Alm Hi	100
SpO ₂ Alm Lo	90
PR Alm Hi	130
PR Alm Lo	50
+/-	+
Exit	

Battery Installation

1. Open the battery compartment by sliding the battery door open.
2. Place the two AAA batteries into the battery compartment ensuring correct position.
3. Slide the battery door back in position ensuring battery compartment is fully closed.

Important:

- To prevent damage ensure correct polarity when installing/replacing batteries.
- To prevent damage install negative (-) side of battery first.
- For extended periods of non-use remove the batteries from the oximeter.

Lanyard Installation

1. Thread the thin end of the lanyard through the hole at the rear of the case.
2. Loop the wider end of the lanyard through the thin end before pulling it tight.

Operation Instructions

Pic 1.

1. Install two AAA batteries (See battery installation).
2. Insert finger into the oximeter as shown in Pic 1. The finger nail must be facing upwards.
3. Switch the device on by pressing the power button.
4. Avoid excessive body movement

during measurement.

5. Measurement data will be displayed on the device screen.
6. For best performance and accurate measuring, apply the HbO-3000 to the finger that provides the highest PI value.
7. Remove finger to stop measurement. Device will automatically switch off after 8 seconds.

Note: Use sterile wipes to clean the inside and outer case of the HbO-3000 oximeter between patients. The HbO-3000 oximeter is constructed of latex free materials.

Product Specifications

1. **Display Type:** OLED (Organic Light Emitting Diode).
2. **SpO₂:**
Measurement range: 35%-100%
Accuracy: 70% - 99% ±3% ≤69%, no

definition

3. Pulse Rate and PI:

Measurement range: 30 - 250 BPM

Accuracy: 30 - 250 BPM $\pm 2\%$

Pulse rate display model: Bargraph
Indicator

PI parameter range: 0-30%. 20%-30%
Unspecified

PI Accuracy: 0.1%

4. Power Requirements:

Two AAA alkaline batteries

Power consumption: Less than 30mA

Low power indication: The low voltage
warning will illuminate when the
voltage is less than 2.6V($\pm 0.2V$)

Battery Life: Two AAA 1.5V, 600mAh
alkaline batteries 20 hours operation

5. Dimensions:

LxWxH: 57mm x 31.5mm x 30.5mm

6. Environment Requirements:

Operating Temperature: 5 - 40°C

Storage Temperature: -20- 70°C

Ambient Humidity:

Operating:15%-90%;

Storage:15%-93%

7. EMC of this product complies with IEC60601-1-2 standards.
8. Measurement performance in low perfusion condition: Required test equipment (FLUKE-INDEX 2XL) the pulse wave is functional without failure when the simulation pulse wave amplitude is set at 0.4%.
9. Interference Resistance Capacity against ambient light: Device works normally. Mixed noise produced by FLUKE-INDEX 2XL.

Product Accessories

1. Lanyard
2. 2 x AAA batteries
3. User manual

Maintenance and Storage

- Replace the 2 x AAA batteries when low voltage warning is indicated.
- Clean the internal/external surfaces of the HbO-3000 fingertip oximeter before and after patient use.
- For extended periods of non-use remove the batteries from the oximeter.
- Storage conditions: Ambient temperature 10°C-40°C(14-104°F) and relative humidity 10%-80%.
- Store the unit in a damp free, dry environment.
- Please comply with local authority regulations when disposing/recycling of the device and batteries.

Calibration Method

Factory calibration is conducted using a FLUKE-INDEX 2XL oxygen analogue meter, and BCI oximeter calibration meter.

Note: The test methods used to

establish the SpO₂ accuracy requires clinical testing. The oximeter is used to measure the arterial hemoglobin oxygen saturation levels and these levels are to be compared to the levels determined from arterial blood sampling with a CO-oximeter.

Troubleshooting

Problems	Possible reason	Solution
SpO ₂ or Pulse Rate not being displayed	Finger is not inserted correctly	Re-insert or try a different finger
SpO ₂ or Pulse Rate not being displayed	Excessive body movement	Reduce and stabilise body movement

<p>The oximeter will not switch on</p>	<p>Exhausted batteries Batteries incorrectly installed</p> <p>Faulty oximeter</p>	<p>Replace batteries Reinstall the batteries</p> <p>Please contact your local customer service centre</p>
<p>Device display turns off suddenly</p>	<p>The product will automatically power off after 8 seconds if no signal is detected</p> <p>Exhausted batteries</p>	<p>Normal</p> <p>Replace the batteries</p>

Symbol Definitions

	Type B applied part
	Attention, refer to the manual before using
%SpO ₂	Oxygen saturation
PR _{bpm}	Heart rate (BPM)
	Low power indicator
 SpO ₂	No SpO ₂ alarm
SN	Serial number

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DeVilbiss Healthcare GmbH
Manufacturer: **Beijing Safe Heart Technology Ltd.**
Shanghai International Holding Corp.
GmbH (Europe)